

TIDEWALL GUNRIG

M	T	SV	W	LD	OC
4"	8	3+	14	7+	0

AUFSTELLUNG

1x Tidewall Gunrig

Equipped with: Supremacy railgun

RANGED WEAPONS

WEAPON	RANGE	A	BS	S	AP	D
Supremacy railgun DEVASTATING WOUNDS, TWIN-LINKED	72"	1	5+	20	-5	D6+6

ABILITIES

RULES - Deadly Demise D3, Firing Deck 11

Fortification

While an enemy unit is only within Engagement Range of one or more Fortifications from your army:

- That unit can still be selected as the target of ranged attacks, but each time such an attack is made, unless it is made with a Pistol, subtract 1 from the Hit roll.
- Models in that unit do not need to take Desperate Escape tests due to Falling Back while Battle-shocked, except for those that will move over enemy models when doing so.

Tidewall Cover

Each time a ranged attack is allocated to a model, if that model is not fully visible to every model in the attacking unit because of this FORTIFICATION, that model has the Benefit of Cover against that attack.

Tidewall Gunrig (Transport)

Capacity: This model has a transport capacity of 11 T'AU EMPIRE INFANTRY models. It cannot transport BATTLESUIT, KROOT or VESPID STINGWINGS models.



WEAPON RULES

Devastating Wounds

Each time an attack made with a [DEVASTATING WOUNDS] weapon results in a critical wound, the attack sequence for that attack ends and the target unit suffers a number of mortal wounds equal to the D characteristic of that weapon. These are inflicted after resolving any normal damage inflicted by those attacks.

Mortal wounds inflicted by [DEVASTATING WOUNDS] weapons can damage a maximum of one model for each critical wound; any remaining mortal wounds inflicted by that attack are lost.

Example: An attack made with a [DEVASTATING WOUNDS] weapon with a D characteristic of 3 results in a critical wound against an Intercessor Squad, so inflicts 3 mortal wounds. The first 2 mortal wounds are sufficient to destroy 1 Intercessor model, so the remaining mortal wound is lost.

Twin-linked

Each time an attack is made with a [TWIN-LINKED] weapon, you can re-roll the wound roll.

